KULLANIM AMACI

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INTENDED USE

(English)

DIACHECK CELL Blood Grouping internal quality control reagents are used to evaluate the performance of tests, systems and users used in blood group testing laboratories and blood centers.

SUMMARY AND PRINCIPLES

It can be used to test all kinds of blood group methods (Gel, Microplate, Tube, Tile).

The reagents are of human origin.

It is supplied ready in appropriate preservatives and at appropriate dilution rates to control the tests routinely used in blood centers.

CONTROL MATERIAL

This whole blood control material contains all of the following:

- Patient-1 Erythrocyte Suspension (1%)
- Donor-1 Erythrocyte Suspension (1%)
- Patient-2 Erythrocyte Suspension (1%)
- Donor-2 Erythrocyte Suspension (1%)
- Patient-1 Serum
- Patient-2 Serum

PRECAUTIONS

1. For IVD use only.
2. Although the serology tests prescribed by WHO were performed and negative samples were used, it should not be forgotten that it is a blood material and may contain infectious virus.
3. It should be treated as medical waste and disposed of.

STORAGE AND STABILITY

It should be stored at 2-8 °C. Stable from unopened to expiration date, 30 days after capping.

PRODUCT DEFAULT

It may be an indication that it is outside the expected values. However, excessive darkening in the color of the product or formation of clots indicates that the product has deteriorated.

USAGE

1. Remove the control vial at 2-8 °C and leave at room temperature for 15 minutes.
2. You can mix it in two different ways.
 - a. Hold the bottle vertically and rotate each bottle between your palms for 15-20 seconds.

- b. Holding the vial at the ends between thumb and finger, continue mixing by rapidly inverting the vial 20

IMMUNOHEMATOLOGY CONTROL

times from end to end using a very quick whirling motion of the wrist.

- c. Analyze immediately after mixing. Device analysis can be performed by inverting the vial 5 times.

3. To analyze control materials:

- a. Perform ABO and D identification tests using Patient-1 and Patient-2 erythrocyte and serum samples.
- b. Perform the CrossMatch tests using Patient-1 with Donor-1 and Patient-2 with Donor-2, Patient-1 Plasma and Patient-2 Plasma samples.
- c. Perform direct Coombs tests using Patient-1 and Patient-2 erythrocyte suspensions.
- d. Perform indirect Coombs tests using Patient-1 and Patient-2 serum samples.
- e. Determine the Rh subgroup using Patient-1 and Patient-2 erythrocyte suspensions.

4. After use, the bottle should be closed tightly and immediately put in the refrigerator.

EXPECTED RESULTS

TEST	PATIENT-1	PATIENT-2
A	Positive	Positive
B	Positive	Negative
AB	Positive	Positive
D	Positive	Positive
A1 CELL	Negative	Negative
B CELL	Negative	Positive
CROSS-MATCH	Positive	Negative
INDIRECT COOMBS	Positive	Negative
DIREKT COOMBS	Negative	Negative
RH FENOTİP	C+ c- E- e+ K-	C+ c+ E+ e+ K-

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1. Marion E. Reid Christine Lomas-Francis THE BLOOD GROUP ANTIGEN FactsBook 2004, Elsevier Ltd.
Medios Medikal Bilişim Danışmanlık Tic. Ltd. Şti.