Ürün bilgisi

24LDL01-UN	Meditest HDL	4x40 mL 2x20 mL
24LDL01-AU	Meditest HDL	4x40 mL 2x20 mL
24LDL01-AB	Meditest HDL	4x40 mL 2x20 mL
24LDL01-ER	Meditest HDL	4x40 mL 2x20 mL
24LDL01-AR	Meditest HDL	4x40 mL 2x20 mL

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Wavelength: 572 nm
Cuvette: 1 cm light path
Temperature: 37°C

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Referanslar

- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.1988;26:783-790
- Rifai N, Warnick GR, McNamara JR, Belcher JD, Grinstead GF, Expected Values Handbook of Laboratory Medicine, Li-hua Zhu 1998
- Frantz Jr ID. Measurement of Low-Density-Lipoprotein Cholesterol in Serum: a Status
- Report. Clin Chem 1992; 38:150-160.
- Kannel, W.B., Castelli W.P., Gordon, T., Cholesterol in the Prediction of Artherosclerotic
- Disease; New Perspectives Based on the Framingham Study, Am. Intern. Med., 90:85 (1979).
- Bachoric P. Measurement of Low-Density-Lipoprotein. 245-263. In: Handbook of
- Lipoprotein Testing (eds. Rifai, Warnick and Dominiczak), 2nd edition, AACC press, 2000.
- National Cholesterol Education Program. Expert Panel on Detection, Evaluation, and
- Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II). NIH Publication No. 93-3095, 1995.
- Naito HK, Strong JP, Scott MG, Roheim PS, Asztalos BF, Zilversmit DB, Srinivasan SR,
- Berenson GS, Wilson PWF, Scanu AM, Malikow MR. Atherogenesis: current topics on etiology and risk factors. Clin. Chem 1995; 41:132-133 No. 1.
- National Committee for Clinical Laboratory Standards, National Evaluation Protocols for
- Interference Testing, Evaluation Protocol Number 7, Vol. 4, No.8, June 1984.
- Wieland H, Seidel D. Quantitative Lipoprotein Electrophoresis.
- Westgard, J.O., Carey, R.N., Wold, S., Criteria for judging precision and accuracy in method development and evaluation.

ClinicalChemistry1974;20:825-833.

- Armstrong V, Seidel D. Evaluation of a Commercial Kit for the Determination of LDL Cholesterol in Serum Based on Precipitation of LDL with Dextran Sulfate. Ärztl. Lab. 1985;31:325-330.
- National Institutes on Health Publication no. 93-3095, September 1993.
- Stone NJ, Robinson JG, Lichtenstein AH, et al. 2013 ACC/AHA guideline on the treatment of blood cholesterol to reduce atherosclerotic cardiovascular risk in adults: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines. J Am Coll Cardiol 2014;63:2889-2934.

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Medios Medikal Bilişim Danışmanlık Tic. Ltd. Şti.